## stratec

## **CONTENT LOADING**

# INDUSTRIALIZATION OF FUNCTIONALIZED CONSUMABLES



spotting of content

### **Application**

many options for customized content loading

#### CONTENT FILLING

- customized content filling or pre-loading of lab-on-a-chip consumables
  - · liquids and solids, planar or localized on chip
  - · variable volumes (pl to ml)
  - · contactless dispensing/spotting
  - · diverse reagents (biological, chemical)
- in combination with suitable packaging
  - · packaging of reagents-kits
  - sterile dust-free shock proof non degassing – waterfree
  - · blister packaging vacuum inert gas filled e.g. N
  - · chilled products and frozen cargo

#### YOUR BENEFITS

#### industrialization

 robust and reproducibility processes for mass production

#### customized solutions

- · design for manufacturing
- prototype test spotting to high volume production

#### optimizing processes

- · large choice of systems (pl to ml)
- content gets optimized on substrate properties
- · high expertise in materials and applications

#### distribution (e.g. frozen cargo)



hybcell consumable for multiplex IVD assays



test-kit for Cube Dx including hycells and reagents